



PROTOCOL 1.

DNA isolation from herbarium leaf tissue using DNeasy Plant Mini Kit[®] and QIAcube[®] robotic workstation

Things to be done before the start:

- All centrifugation steps are carried out at room temperature (15–25°C)
- Buffer AP1 and Buffer AW1 concentrate may form precipitates upon storage. If necessary, warm to 65°C to redissolve (before adding ethanol to Buffer AW1). Do not heat Buffer AW1 after ethanol has been added
- Before using for the first time, add the appropriate amount of ethanol (96–100%) to Buffers AW1 and AW2 as indicated on the bottles to obtain working solutions.
- The QIAcube performs fully automated processing of up to 12 samples.
- Pour AW1, AW2 and AE buffers into three 30 ml QIAcube reagent bottles.
- Place 30 ml QIAcube reagent bottles into the position 1 (AW1), 2 (AW2) and 3 (AE) of the QIAcube reagent bottles rack.
- Place DNeasy Mini spin column (white), QIAshredder spin column (violet column with cut off lid) from DNeasy[®] Plant Mini kit and 1,5 ml microcentrifuge tube into the position 1, 2 and 3 of the QIAcube[®] Rotor Adapter.
- Place the loaded QIAcube Rotor Adapter into the QIAcube centrifuge bucket.
- Place two prefilled tip racks onto the worktable of the QIAcube.

- 1) Place ca. 20 mg of herbarium leaf tissues into a 2 ml safe lock (Eppendorf[®]) microcentrifuge tubes. Place the tubes into the TissueLyser Adapter Set 2 x 24 (Qiagen[®]). Add 5 mm autoclaved stainless steel bead into each microcentrifuge tube.
- 2) Fix TissueLyser Adapter Set 2 x 24 into the clamps of the TissueLyser (Qiagen[®]) and grind the samples for 1 min at 30 Hz.
- 3) Add 400 µl Buffer AP1, 4 µl 2-mercaptoethanol and 4 µl RNaseA stock solution into each microcentrifuge tube with disrupted leaf tissue and vortex vigorously using vortex mixer.
 - Carry out this step in a fume hood
- 4) Incubate the mixtures in microcentrifuge tubes for 10 min at 65°C. Mix 2-3 times during incubation using vortex mixer.
 - This step lyses the cells
- 5) Add 130 µl Buffer P3 to each lysate, mix, and incubate for 7 min on ice
 - This step precipitates detergent, proteins, and polysaccharides
- 6) Centrifuge the lysates for 5 min at 20,000 x g
- 7) Transfer the supernatants to a new 2 ml safe lock (Eppendorf[®]) microcentrifuge tubes.
 - Be careful not to pick up pellet at the bottom of the tube during pipetting process
- 8) Open the QIAcube door and place 2 ml tubes into the QIAcube shaker adapter. Place tube lids into the slot at the edge of the shaker adapter.
 - Remove caps from 30 ml QIAcube reagent bottles



- 9) Close the QIAcube door and launch QIAcube robotic workstation choosing protocol: Purification of total DNA from plant cells and tissues (see QIAcube Protocol Sheet).
 - In these steps DNA molecules are bounded, washed, dried and finally eluted from DNeasy Mini spin column as a pure DNA isolate.
- 10) When the protocol run has finished, a message is displayed in the touchscreen confirming that the samples have been processed.
- 11) Remove the 1,5 ml microcentrifuge tubes containing purified DNAs from the rotor adapters.
 - Discard sample tubes, used rotor adapters, and reagent according to safety regulations.
 - Replace the lids of the reagent bottles and close tightly
 - Empty the waste drawer
- 12) Isolated DNA solution should be kept at -20 °C for longer period or at 4-8 ° C for several days

References:

- 1) DNeasy® Plant Handbook
- 2) QIAcube® User Manual
- 3) QIAcube and DNeasy Plant Mini Kit - QIAcube® Protocol Sheet



Laboratory equipment and accessories used for DNA isolation



TissueLyser (Qiagen®)



Dneasy® Plant mini kit



GVLab vortex mixer (Gilon®)



Bio TDB-100 Dry Block Biosan®



Centrifuge Hettich® Mikro 200



QIAcube® robotic workstation



QIAcube® Filter-Tips 1000 µl



QIAcube® 30 ml Reagent Bottles



QIAcube® Rotor Adapter

Plant material

