



PROTOCOL 2.

DNA isolation from 124 yers old Balkan chamois bone remains using QIAquick[®] PCR purification kit and QIAcube[®] robotic workstation

Things to be done before the start:

- All surfaces used in the laboratory should be wiped down with a 10% of sodium hypochlorite and 70% ethanol solution before and after each use.
- Equipment (e.g. stainless steel grinding jars and balls, pipettes, tips) should be irradiated under the UV light for 30 minutes, before and after use
- Add the appropriate amount of ethanol (96–100%) to Buffer PE before use to obtain working solutions.
- The QIAcube performs fully automated processing of up to 12 samples.
- Pour PE and EB buffers into two 30 ml QIAcube reagent bottles.
- Place 30 ml QIAcube reagent bottles into the position 5 (PE) and 6 (EB) of the QIAcube reagent bottles rack.
- Place QIAquick column and 1,5 ml microcentrifuge tubes into the position 1 and 3 of the QIAcube[®] Rotor Adapters.
- Place the loaded QIAcube Rotor Adapters into the QIAcube centrifuge buckets.
- Place two prefilled tip racks (1000 µl) onto the worktable of the QIAcube.

- 1) Cut 0.5 to 1 g of molar and/or maxilla bone using Dremel[®] machine.
- 2) Wash bone cut twice with 10% of sodium hypochlorite for 30 s, soak in water for 5 min, rinse using 96% ethanol and then dry at room temperature.
- 3) Ground bone cut for 20 min at 30 Hz in the TissueLyser (Qiagen[®]) using two stainless steel grinding jars and 10 mm stainless steel grinding balls.
- 4) Incubate 0,5 g of bone powder in 10 ml of lysis buffer (0,5 M EDTA pH 8, 1 % N-laurylsarcosinate and 6,7 mg of proteinase K) overnight in an incubator at 56 °C on a rocking or rotated platform.
- 5) Centrifuge lysates with completely digested bone powder for 5 min at 1800 x g to pellet any particulates.
- 6) Concentrate lysates to approximately 300 µl by centrifugation for 50 to 60 min at 5000 x g in Amicon[®] Ultra-15 100K Centrifugal Filter Devices.
- 7) Transfer each 300 µl of concentrated lysates to a new 2 ml safe lock (Eppendorf[®]) microcentrifuge tubes.
- 8) Add 1500 µl Buffer PB to each concentrated lysate and mix.
- 9) Open the QIAcube door and place 2 ml tubes into the QIAcube shaker adapter. Place tube lids into the slot at the edge of the shaker adapter.
 - Remove caps from 30 ml QIAcube reagent bottles
- 10) Close the QIAcube door and launch QIAcube robotic workstation choosing: QIAcube[®] Customized Protocol Sheet IDC1179 (see QIAcube Protocol Sheet).



- In these steps DNA molecules are bounded, washed, dried and finally eluted from DNeasy Mini spin column as a pure DNA isolate.
- 10) When the protocol run has finished, a message is displayed in the touchscreen confirming that the samples have been processed.
 - 11) Remove the 1,5 ml microcentrifuge tubes containing purified DNAs from the rotor adapters.
 - Discard sample tubes, used rotor adapters, and reagent according to safety regulations.
 - Replace the lids of the reagent bottles and close tightly
 - Empty the waste drawer
 - 12) Isolated DNA solution should be kept at -20 °C for longer period of time or at 4-8 ° C for several days

References:

- 1) Amicon® Ultra15 Centrifugal Filter Unit - User guid
- 2) Amory et al. 2012. Automatable full demineralization DNA extraction procedure from degraded skeletal remains. Forensic Science International: Genetics 6: 398–406.
- 3) QIAcube® - User Manual
- 4) QIAcube and QIAquick PCR Purification Kit - QIAcube® Protocol Sheet
- 5) QIAquick® - PCR Purification Kit



6) Laboratory equipment and accessories used for DNA isolation



Dremel 3000-15 Multitool



Tissue Lyser Qiagen®



Grinding jars and balls
for Tissue Lyser



Stuart® Hybridisation
oven SI30H



Amicon® Ultra 15 mL
centrifugal filter devices



Sigma 4-16 Centrifuges



QIAquick® PCR purification kit



QIACube®



QIACube® Filter-Tips 1000 µl



QIACube® 30 ml
Reagent Bottles



QIACube® Rotor Adapter